

## Bulk pigment extraction from lake sediments

### Consumables and equipment

#### Chemical reagents

- Acetone **100%** (grade **HPLC**): solvent.
- Acetone **technical**: cleaning.
- **Redistilled water** (*MilliQ*).
- Technical nitrogen.
- Neodisher Laboclean.

Neodisher LaboClean FLA: liquid, highly-alkaline cleaning, dispergator agent.

#### Equipment

- Laboratory fume hood.

#### Acetone

Acetone ( $C_3H_6O$ ) is an organic solvent requiring work under the **running fume hood**. Acetone is highly **flammable** and is an **irritant**.

- Ultrasound bath.
- Centrifuge.
- Vortex shaker.
- Laboratory furnace for glassware baking.
- Laboratory oven.
- Freeze drier.
- Glass beakers for acetone pipetting.
- Acetone pipette **1000  $\mu$ l** and tips (for acetone).
- Heating plate.
- Nitrogen evaporator.

- Needles for N<sub>2</sub> evaporator.
- Metal rack for amber vials on the heating plate.

For example *Turbovap* instrument.

- Racks for centrifuge tubes and amber vials.
  - Analytical balance, antistatic device and stainless steel spatula.
  - Ice bath.
- Plastic tray with cooling packs.
- Crimper (for closing HPLC vials).

## Other

- Polypropylene (PP) centrifuge tubes (falcons) **25 ml** or **50 ml** (Corning™).
- Other plastics will melt in contact with acetone.  
One-time use.
- Amber vials with caps **25 ml** or **50 ml**.
- Tape the labels.
- PTFE hydrophobic filters of **0.2 μm** (**13 mm**).
- One-time use.
- Polypropylene (PP) small syringes with Luer Lock for filtering (**2 ml**).
- One-time use.

### ! Important

Always check size and fit compatibility.

- Aluminum foil.
  - Glass Pasteur pipettes.
- One-time use.
- Silicone nipples for pipettes.

## Cleaning

### Amber vials

- Put into soap water with Neodisher Laboclean for few hours.
- Wash in the washing machine.  
Glassware program and Neodisher Laboclean.
- Dry in the drying oven.
- Close with aluminum foil (cleaned with technical acetone).
- Bake in the furnace in **500 °C**.

Minimum **3 h**.

### Amber vials caps

- Rinse with tap water.
- Put into beaker filled with **technical acetone** and **redistilled water (MilliQ)** mixture (**1:1**).
  - Sonicate for **10 minutes**.
- Put into beaker filled with **technical acetone**.
  - Sonicate for **10 minutes**.
- Dry in the drying oven.

Washing solvents can be used again. Store in the glass bottles.

### Evaporator needles

- Put into beaker filled with **technical acetone**.
- Sonicate **3 × 10 minutes**.

### Disposal

- Acetone accordingly.
- PP falcons and glass Pasteur pipettes need to be washed and recycled.

## Sediment preparation

- Freeze dry.  
High temperature speeds up pigment degradation.
- Loss on ignition **550 °C**.  
Paleoenvironmental Research Laboratory: calculate organic matter from CNS.
- Weigh **0.5 g** homogenized sediments into polypropylene falcons (PP) **25 ml**.  
Clean tools with technical acetone, use antistatic device.  
If there is not enough material, use no less than **0.25 g**.
- Label falcon and cap with sample name.

## Pigment extraction

- Prepare workplace, cover surfaces under the fume hood with aluminum foil cleaned with technical acetone.
- Switch off the lights whenever possible.
- Fill glass beaker with **100%** acetone.  
Don't pipette from the bottle.
- Rinse pipette with acetone.  
Pipette into waste beaker.
- Add **5 ml** of acetone to the falcon.  
pipette right, there should be no air bubbles in the tip.
- Shake with vortex for **1 minute**.
- Sonicate for **1 minute**.
- Centrifuge:
  - **10 minutes**.
  - **5000 rpm**.

This step needs to be adjusted to the centrifuge. Lower rpm requires longer centrifugation.

- Check if supernatant is clear after centrifugation.  
If it is not clear centrifugate again.
- Remove (transfer) supernatant using glass Pasteur pipette into amber vial.  
Label amber vial accordingly.  
Be sure all supernatant is transferred.
- Repeat two times using **5 ml** acetone each time.
- Hold labelled Pasteur pipettes in a glass beaker between the extraction steps.  
Every amber vial and falcon has its own Pasteur pipette.

- Amber vials between extraction steps should be covered from light and stored in the ice bath.

After **three extraction cycles**, if supernatant still has a strong color continue with **additional extraction steps**.

It is important that **all samples** are prepared **the same way**.

Extracts can be stored **for a few days** in the **-20 °C** temperature.

## Concentration and reconstitution

### Evaporation

- Set heating plate and nitrogen evaporator.
- Evaporate samples to dry residue.

Before use be sure that needles are clean (Section ??).

- Put amber vials in the rack on the heating plate set to **35 °C**.
- Open nitrogen flow.

Set flow so small flexion of the liquid surface is visible.

- Evaporate to dry residue (around **4 to 6 h**).

### Reconstitution and filtering

- Add **2 mL 100%** acetone to the amber vial using **1000 µl** pipette.

Be very precise, check for air bubbles in the pipette tip.

- Homogenize solution with the glass Pasteur pipette.
- Rinse the vial walls and *vortex*.

All pigments need to redissolve in the acetone.

- Take a **2 mL** syringe, attach the **0.2 µm** PTFE hydrophobic filter on.
- Pipette directly with a Pasteur pipette into the syringe.
- Pass the solvent through the filter into the small HPLC glass vials with the piston.
- Loose the filter from the syringe to get all solvent out.
- Close the HPLC vials with a crimper and store your samples in a labelled sample box at -20 °C.

## Credits

Andrea Sanchini, Giulia Wienhues and Paul Zander – University of Bern, Switzerland.

## Changelog

- 09.12.2022, MZ – initial version.
- 09.02.2023, MZ – full version.
- 17.03.2023, MZ – update with filtering and callouts.

[maurycy.zarczynski@ug.edu.pl](mailto:maurycy.zarczynski@ug.edu.pl)

<http://geomorfologia.ug.edu.pl>