

Bulk pigment extraction from lake sediments

Consumables and equipment

Chemical reagents

- Acetone **100%** (grade **HPLC**): solvent.
- Acetone **technical**: cleaning.
- **Redistilled water** (*MilliQ*) (or substitute).
- Technical nitrogen.
- Neodisher Laboclean (or substitute).
Neodisher LaboClean FLA: liquid, highly-alkaline cleaning, dispersant agent.

Equipment

- Laboratory fume hood.

Acetone

Acetone ($\text{C}_3\text{H}_6\text{O}$) is an organic solvent requiring work under the **running fume hood**. Acetone is highly **flammable** and is an **irritant**.

- Ultrasound bath.
- Centrifuge.
- Vortex shaker.
- Laboratory furnace for glassware baking.
- Laboratory oven.
- Freeze drier.
- Glass beakers for acetone pipetting.
- Acetone pipette **1000 μL** and tips (for acetone).
- Heating plate.
- Nitrogen evaporator.
 - Needles for N_2 evaporator.
 - Metal rack for amber vials on the heating plate.

For example *Turbovap* instrument.

- Racks for centrifuge tubes and amber vials.
- Analytical balance, antistatic device and stainless steel spatula.
- Ice bath.

Plastic tray with cooling packs.

Other

- Polypropylene (PP) centrifuge tubes (falcons) **25 mL** or **50 mL** (Corning™).

Other plastics will melt in contact with acetone.

One-time use.

- Amber vials with caps **5 mL**.

Tape the labels.

Two sets are needed, for unfiltered and filtered solvent.

- Polypropylene (PP) small syringes with Luer Lock for filtering (**5 mL**).

One-time use.

- PTFE hydrophobic filters of **0.2 m**.

One-time use.

! Important

Always check size and fit compatibility for syringes and filters.

- Aluminum foil.
- Glass Pasteur pipettes.

One-time use.

- Silicone nipples for pipettes.

Cleaning

Amber vials

- Put into soap water with Neodisher Laboclean for few hours.
- Wash in the washing machine.

Glassware program and Neodisher Laboclean.

If dishwasher is unavailable do it manual.

- Dry in the drying oven.
- Close with aluminum foil (cleaned with technical acetone).
- Bake in the furnace in **380 °C**.

! Baking

Minimum **3 h**.

Never bake measurement glass.

Amber vials caps

- Rinse with tap water.
- Put into beaker filled with **technical acetone** and **redistilled water (MilliQ)** mixture (**1:1**).
 - Sonicate for **10 minutes**.
- Put into beaker filled with **technical acetone**.
 - Sonicate for **10 minutes**.
- Dry in the drying oven.

Washing solvents can be used again. Store in the glass bottles.

Evaporator needles

- Put into beaker filled with **technical acetone**.
- Sonicate **3 × 10 minutes**.

Disposal

- Acetone accordingly.
- PP falcons and glass Pasteur pipettes need to be washed and recycled.

Sediment preparation

- Freeze dry.
 - High temperature speeds up pigment degradation.
- Loss on ignition **550 °C** or **CNS** to calculate organic matter.
- Weigh **0.5 g** homogenized sediments into polypropylene falcons (PP) **25 mL**.
 - Clean tools with technical acetone, use antistatic device.
 - If there is not enough material, use no less than **0.25 g**.
 - With minimal **OM** or **TOC** concentrations weigh **1.00 g** and more if needed.
- Label falcon and cap with sample name.

Pigment extraction

- Prepare workplace, cover surfaces under the fume hood with aluminum foil cleaned with technical acetone.
- Switch off the lights whenever possible.
- Fill glass beaker with **100%** acetone.
Don't pipette from the bottle.
- Rinse pipette with acetone.
Pipette into waste beaker.
- Add **5 mL** of acetone to the falcon.
pipette right, there should be no air bubbles in the tip.
- Shake with vortex for **1 minute**.
- Sonicate for **1 minute**.
- Centrifuge:
 - **10 minutes**.
 - **5000 rpm**.

This step needs to be adjusted to the centrifuge. Lower rpm requires longer centrifugation.

- Check if supernatant is clear after centrifugation.
If it is not clear centrifugate again.
- Remove (transfer) supernatant using glass Pasteur pipette into amber vial.
Label amber vial accordingly.
Be sure all supernatant is transferred.
- Repeat two times using **5 mL** acetone each time.
- Hold labelled Pasteur pipettes in a glass beaker between the extraction steps.
Every amber vial and falcon has its own Pasteur pipette.
- Amber vials between extraction steps should be covered from light and stored in the ice bath.

After **three extraction cycles**, if supernatant still has a strong color continue with **additional extraction steps**.

It is important that **all samples** are prepared **the same way**.

Extracts can be stored **for a few days** in the **-20 °C** temperature.

Concentration and reconstitution

Evaporation

- Set heating plate and nitrogen evaporator.
- Evaporate samples to dry residue.
Before use be sure that needles are clean.
- Put amber vials in the rack on the heating plate set to **35 °C**.
- Open nitrogen flow.
Set flow so small flexion of the liquid surface is visible.
- Evaporate to dry residue (around **4 to 6 h**).
- Prepare blank sample from acetone.

Reconstitution and filtering

- Add **5 mL 100%** acetone to the amber vial #1 (unfiltered) using **1000 µL** pipette.

Be very precise, check for air bubbles in the pipette tip.

- Homogenize solution with the glass Pasteur pipette.
- Rinse the vial walls and *vortex*.

All pigments need to redissolve in the acetone.

- Take a **5 mL** syringe, attach the **0.2 µm** PTFE hydrophobic filter on.
- Pipette from the amber vial #1 (unfiltered) directly with a Pasteur pipette into the syringe.
- Pass the solvent through the filter into the amber vial #2 (filtered).
- Loose the filter from the syringe to get all solvent out.
- Close the vial with a cap and store at **-20 °C**.

Credits

Andrea Sanchini, Giulia Wienhues and Paul Zander – University of Bern, Switzerland.

Changelog

- 09.12.2022, MZ – initial version.
- 09.02.2023, MZ – full version.
- 17.03.2023, MZ – update with filtering and callouts.
- 27.03.2023, MZ – NAU final version 1.

maurycy.zarczynski@ug.edu.pl